

# AMINOACIDS

Column: **CGB-50ht**, P/N CG-531332  
Size: 30m x 0.25mm x 0.15µm  
Injection: 2 µl standard AA-S-18 Sigma (2.5µmol/ml), split, 300°C  
Carrier gas: He, 1mL/min  
Program temperature: 50°C (2min) @ 30°C/min to 350°C (3min)  
Detector: MS Polaris Q, EI, 200°C, transfer line 200°C

*Chromatogram provided by Antonio Tintó from Moehs S.A.*

## Peak Name

- 1 Alanine
- 2 Glycine
- 3 Valine
- 4 Leucine
- 5 Isoleucine
- 6 Serine
- 7 Threonine
- 8 Proline
- 9 Hydroxyproline
- 10 Methionine
- 11 Phenylalanine
- 12 Lysine
- 13 Histidine
- 14 Tyrosine
- 15 Cystine

