

ANALYSIS OF FIRE RESIDUES (PYROLISIS)

Column: **CGB-5**, P/N CG-120262

Dimensions: 60m x 0.25mm x 0.25 μ m


Injection: 1 μ L, splitless

Carrier gas: He, 1.8 mL/min

Oven temperature: 150°C @ 2°C/min to 225°C

Detector: FID, 300°C

*Chromatogram provided by Montse Elias and Jordi Codina
from Laboratori General d'Assaigs i Investigacions*



Peak Name

- 1- Toluene
- 2- Benzene
- 3- Xylene

