

PUFA III

Column: **CronusWax-280**, 15m x 0.10mm x 0.10 μ m (P/N: CG-830111)


Injection: 280°C, split 200:1, precision liner

Detector: FID, 280°C

Carrier Gas: H₂, 45 psi

Oven: 100°C (0.5min) to 280°C @ 50°C/min (2min)

Sample: 0.1 μ L PUFA III – Partially Hydrogenate Menhaden Oil diluted to 100mg/ml in hexane



Peak Name

1. C14:0
2. C16:0
3. C16:1v 7
4. C16:2v 4
5. C16:3v 4
6. C16:4v 1
7. C18:0
8. C18:1v 9
9. C18:1v 7
10. C18:2v 6
11. C18:3v 4
12. C18:3v 3
13. C18:4v 3
14. C20:1v 9
15. C20:4v 6
16. C20:4v 3
17. C20:5v 3
18. C22:5v 3
19. C22:6v 3