

## SEPARATION OF PESTICIDES

Column: **Cronus X5**, P/N CG-820232

Dimensions: 30m x 0.25mm x 0.25  $\mu$ m

Injection: 1  $\mu$ L standard (600 mg/L), split 1:20, 200°C

Carrier gas: He, 1mL/min

Oven program: 90°C(1min) @ 20°C/min to 200°C @ 3°C/min to 220°C

Detector: ECD, 250°C

*Chromatogram provided by AINIA*

	<u>Peak Name</u>	<u>RT (min)</u>
1	Nonadecane (Internal Standard)	12.79
2	Chlorpiryphos	14.59

