

SEPARATION OF SOLVENTS

Column: **CGB-624**, P/N CG-603075


Dimensions: 75m x 0.53mm x 3.0 μ m

Injection: 0.2 μ L, split 1:5, 260°C

Carrier gas: H₂, constant pressure 7.8 psi (53.74 KPa).

Oven temperature: 40°C(5min) @ 7°C/min to 240°C

Detector: FID, 280°C



Peak Name

- 1- Heptane
- 2- Acetone
- 3- Tetrahydrofurane
- 4- Ethyl acetate
- 5- Isopropyl acetate
- 6- Methanol
- 7- MEK
- 8- Isopropanol
- 9- Ethanol
- 10- n-Propyl acetate
- 11- MIBK
- 12- n-Propanol
- 13- Toluene
- 14- Butyl acetate
- 15- Methoxypropanol
- 16- Methoxypropyl acetate