

# SEPARATION OF SOLVENTS

Column: **CGB-WAX**, P/N CG-141253

Dimensions: 50m x 0.32mm x 1.2  $\mu$ m

Injection: 1  $\mu$ L standard (500 ng/mL comp.), split 1:25, 260°C

Carrier gas: He, constant pressure 12 psi (82.7 Kpa)

Oven temperature: 65°C(7min) @ 4°C/min to 117°C

Detector: FID, 260°C

*Chromatogram provided by  
Jaume Piedrabuena from Danisco*

## Peak Name

- 1 Heptane (isomers mixture)
- 2 Acetone
- 3 Tetrahydrofurane
- 4 Ethyl acetate
- 5 Isopropyl acetate
- 6 Methanol
- 7 MEK
- 8 Isopropanol
- 9 Ethanol
- 10 n-Propyl acetate
- 11 MIBK
- 12 n-Propanol
- 13 Toluene
- 14 Butyl acetate
- 15 Methoxypropanol
- 16 Methoxypropyl acetate

