

MIXTURE OF SOLVENTS AND ISOMERS OF N-HEPTANE

Column: **CGB-624**, P/N CG-603075


Dimensions: 75m x 0.53mm x 3.0 μ m

Injection: mixture of solvents (wet needle), split 1:100, 250°C

Carrier gas: He, constant pressure 8 psi (55.7 Kpa)

Oven temperature: 40°C(15 min) @ 15°C/min to 75°C(15 min)

Detector: FID, 250°C



Peak Name

- 1 Methanol
- 2 Acetone
- 3 Isopropanol
- 4 Ethyl acetate
- 5 Tetrahydrofurane
- 6 Isomers n-heptane