

REACTION PRODUCTS OF 1,1-DIETHOXY BUTANE

Column: **CGB-624**, P/N CG-603035

Size: 30m x 0.53mm x 3.0 μ m

Injection: 0.5 μ l, split 1:5, 260 $^{\circ}$ C

Carrier Gas: He, 6 psi

Program temperature: 40 $^{\circ}$ C (6min) @ 30 $^{\circ}$ C/min to 200 $^{\circ}$ C (5min)

Detector: FID, 260 $^{\circ}$ C

Chromatogram provided by Ion Aguirre from Escuela Superior de Ingeniería de Bilbao (Spain)

Peak Name

- 1- Butanal
- 2- Ethanol
- 3- 1,1-Diethoxy butane
- 4- THF

