

## FREE ACIDS IN WATER

Column: **CGB-FFAP**, P/N CG-151035

Dimensions: 30m x 0.53mm x 1.0  $\mu$ m

Injection: 1  $\mu$ L, split

Carrier gas: He, 4 psi (27.56 KPa)

Oven temperature: 120°C @ 4°C/min to 220°C

Detector: FID, 275°C

### Peak Name

- 1- Acetic acid
- 2- Propionic acid
- 3- Isobutyric acid
- 4- Butyric acid
- 5- Isovaleric acid
- 6- Valeric acid
- 7- Isocaproic acid
- 8- Caproic acid
- 9- Heptanoic acid

