

## HYDROCARBONS

Column: **CGB-1ht** P/N CG-610133

Size: 30m x 0.32mm x 0.1µm

Injection: hydrocarbons standard 1250 ppb, splitless, 250°C

Carrier gas: He, constant flow 2 mL/min

Program Temperature: 50°C(1 min) @ 40°C/min to 320°C(5 min)

Detector: MS, ion source 200°C, Interfase 280°C, scan 20-600

Peak Name

1	C10
2	C12
3	C24
4	C28
5	C30
6	C40

*Chromatogram provided by Vanesa Riu de ILERSAP, Mollerussa (Lleida).*

