

PESTICIDES IN WINES

Column: **CGB-5**, P/N CG-120232

Dimensions: 30m x 0.25mm x 0.25 μ m

Injection: 1 μ L standard in Cyclohexane, 250°C

Carrier gas: He, 1 mL/min

Oven temperature: 145°C @ 0.55°C/min to 220°C

Detector: ECD, 300°C

*Chromatogram provided by J. García from INCAVI,
Vilafranca del Penedès (Barcelona)*

Peak Name

- 1- Methyl parathion
- 2- Vinclozoline
- 3- Phenitrothion
- 4- Diclofluanide
- 5- Chlorpiryphos
- 6- Captan
- 7- Folpet
- 8- Chlozolate
- 9- Procimidone
- 10- Captafol
- 11- Iprodione
- 12- Bromopropilate
- 13- Phenarimol

