

INDUSTRIAL SOLVENTS

Column: **CGB-WAX**, P/N CG-141233

Dimensions: 30m x 0.32mm x 1.2 μ m

Injection: 0.1 μ L, split

Carrier gas: He, 12 psi (82.7 KPa)

Oven temperature: 40°C @ 1°C/min to 70°C @ 7.5°C/min to 125°C

Detector: FID, 250°C

Peak Name

- 1- Acetone
- 2- Methanol
- 3- Ethanol
- 4- p-Xylene
- 5- m-Xylene
- 6- o-Xylene

