

# USP SOLVENTS <USP> COLUMN TRB-G27+ GUARD COLUMN 5M

Column: **CGB-G27**, P/N CG-175035

Dimensions: 30m x 0.53mm x 5.0  $\mu$ m

Oven temp.: 35°C(5') to 175°C@ 8°C/min. to 260°C (16')@35°C/min.

Carrier gas: He, 4.5 psi (31 KPa), 35 cms. to 35°C

Injector temp: 70°C

FID temp: 260°C

Injection: Direct injection of 1  $\mu$ l (Uniliner), standard dissolution in distilled water (1:10)

	<u>Standard</u>	<u>Concentration</u>
1-	Methylene chloride	600ppm
2-	Chloroform	60ppm
3-	Benzene	2ppm
4-	Trichloroethylene	80ppm
5-	1,4 - Dioxan	380ppm

