

## SOLVENTS IN WATER (100 ppm)

Column: **CGB-WAX**, P/N CG-142065

Dimensions: 60m x 0.53mm x 2.0  $\mu$ m

Injection: 1  $\mu$ L, split

Carrier gas: He, 14 psi (96.46 KPa)

Oven temperature: 60°C (Isothermal)

Detector: FID, 280°C

### Peak Name

- 1- Acetone
- 2- Ethyl acetate
- 3- Methanol
- 4- Isopropanol
- 5- Ethanol
- 6- Chloroform

