

SOLVENTS IN WATER

Column: **CGB-5**, P/N CG-120232

Dimensions: 30m x 0.25mm x 0.25 m

Injection: 1 μ L, split

Carrier gas: He

Oven temperature:

Detector: FID

*Chromatogram provided by J. Teixidor and E. Bosch from
Laboratory Dr. Riera*

Peak Name

- 1- Butyl acetate
- 2- Ethyl benzene
- 3- m,p-Xylene
- 4- o-Xylene
- 5- Trimethyl benzenes
- 6- n-Decane
- 7- Tetramethyl benzenes

