

SEPARATION OF SUGARS (AS ALDITOL ACETATES)

Column: **CGB-225**, P/N CG-250232

Dimensions: 30m x 0.25mm x 0.25 μ m

Injection: 1 μ L Sugars standard, split

Carrier gas: H₂, 11 psi (75.8 KPa)

Oven temperature: 220°C (Isothermal)

Detector: FID, 250°C

Peak Name

- 1- Rhamnitol
- 2- Fucitol
- 3- Ribitol
- 4- Arabinitol

