


## Quality Lab Accessories L.L.C

100 Emlen Way #108, Telford, PA. 18969  
Tel. (908) 685-7500 Fax. (908) 526-4151

# *Certificate of Conformance*

<b>Part Number</b>	BSK040-01		
<b>Part Description</b>	#40 Mesh Dissolution Basket, Agilent/Vankel/Varian		
<b>Serial Numbers</b>			

Quality Lab Accessories L.L.C Certifies, Warrants and Represents that the items described herein conform to the specifications set forth in U.S. Pharmacopoeia, General Chapter, (711) Dissolution as last revised, as well as ASTM E 2503-07 and that all components of the listed part number this represents were designed, manufactured or fabricated within strict compliance of these same standards.



No.	USP	QLA
<b>1</b>	20.21 to 20.29 mm	
<b>2</b>	36.52 to 37.03 mm	
<b>3</b>	26.00 to 28.00 mm	
<b>4</b>	21.71 to 22.73 mm	
<b>5</b>	20.21 to 20.29 mm	
<b>6</b>	25.14 to 25.65 mm	
<b>A</b>	Max run-out +/- 1.0 mm	
<b>B</b>	0.25 mm wire diameter with wire openings of 0.38 mm	

Material
<b>316L Stainless Steel</b>

\_\_\_\_\_  
Certified By

\_\_\_\_\_  
Signature

\_\_\_\_\_  
Date