


# Basket Sinker Specifications

## *VanKel Style*



### **Internal Size:**

**Inside diameter:**----- 20.2 mm (0.80 in)  
**Inside height:**----- 11.7 mm (0.46 in)

### **Material Specification:**

**Frame:**----- 316 stainless steel  
**Retainer:** ----- Fluoropolymer (PTFE)  
**Finish:** ----- Electropolish

### **Mesh Specification:**

Catalog Number	VanKel P/N	Mesh	Wire Size (inches)
SNK008-VK	12-3062	8 Mesh, 316 SS, wire cloth	.028 dia
SNK010-VK	12-3060	10 Mesh, 316 SS, wire cloth	.025 dia
SNK040-VK	12-3064	40 Mesh, 316 SS, wire cloth	.010 dia
SNKSPC-VK	custom	SPECIFY MESH SIZE	