


Capsule Sinker Specifications Helix Coil Type



CAPWHT-02

Material	Music Wire
Finish	PTFE Coating
Inside Length	.840 in
Inside Diameter	.370 in
Wire Diameter	.050 in
Pitch	.156 in
End Coil Inside Dia.	.175 in
Maximum Width	.64 in
No. of Coils	8.5
Weight	2.05 gm
Maximum Capsule Size	0

Ordering Information:

QLA Capsule Sinks come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.