


## *Capsule Sinker Specifications Helix Coil Type*



### **CAPWHT-04**

Material	Music Wire
Finish	PTFE Coating
Inside Length	.840 in
Inside Diameter	.385 in
Wire Diameter	.043 in
Pitch	.280 in
End Coil Inside Dia.	.175 in
Maximum Width	.60 in
No. of Coils	5
Weight	1.70 gm
Maximum Capsule Size	0

#### **Ordering Information:**

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.