


Capsule Sinker Specifications

Helix Coil Type



CAPWHT-2S

| | |
|----------------------|---------------------|
| Material | 316 Stainless Steel |
| Finish | Electropolish |
| Inside Length | .840 in |
| Inside Diameter | .370 in |
| Wire Diameter | .043 in |
| Pitch | .156 in |
| End Coil Inside Dia. | .190 in |
| Maximum Width | .64 in |
| No. of Coils | 8.5 |
| Weight | 2.85 gm |
| Maximum Capsule Size | 0 |

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.