


Capsule Sinker Specifications Helix Coil Type



CAPWHT-4S

Material	316 Stainless Steel
Finish	Electropolish
Inside Length	.840 in
Inside Diameter	.385 in
Wire Diameter	.047 in
Pitch	.280 in
End Coil Inside Dia.	.190 in
Maximum Width	.60 in
No. of Coils	5
Weight	1.95 gm
Maximum Capsule Size	0

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.