


Capsule Sinker Specifications

Helix Coil Type



CAPWHT-XL

Material	Music Wire
Finish	PTFE Coating
Inside Length	1.150 in
Inside Diameter	.465 in
Wire Diameter	.043 in
Pitch	.320 in
End Coil Inside Dia.	.175 in
Maximum Width	.69 in
No. of Coils	5.50
Weight	2.15 gm
Maximum Capsule Size	000

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.