


## *Capsule Sinker Specifications*

### *Sotax (Manual) Style*



### **CAPWST-15**

Material	316 Stainless Steel
Finish	Electropolished
Inside Length	15.5 mm
Inside Width	5.0 mm
Wire Diameter	.42 mm
Coil F/L	16.2 mm
Coils	6
Weight	.60 gm
Maximum Capsule Size	5
Sotax Part Number	8280

#### **Ordering Information:**

QLA Capsule Sinker come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.