


Capsule Sinker Specifications

Sotax (Manual) Style



CAPWST-18

Material	316 Stainless Steel
Finish	Electropolished
Inside Length	18.0 mm
Inside Width	6.0 mm
Wire Diameter	.42 mm
Coil F/L	18.8 mm
Coils	8.7
Weight	.75 gm
Maximum Capsule Size	3
Sotax Part Number	7346

Ordering Information:

QLA Capsule Sinker come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.