


## *Capsule Sinker Specifications*

### *Sotax (Manual) Style*



#### **CAPWST-19**

Material	316 Stainless Steel
Finish	Electropolished
Inside Length	19.3 mm
Inside Width	7.0 mm
Wire Diameter	.42 mm
Coil F/L	19.7 mm
Coils	8.0
Weight	1.0 gm
Maximum Capsule Size	2
Sotax Part Number	8283

#### **Ordering Information:**

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.