


Capsule Sinker Specifications

Sotax (Manual) Style



CAPWST-23

Material	316 Stainless Steel
Finish	Electropolished
Inside Length	23.0 mm
Inside Width	8.0 mm
Wire Diameter	.42 mm
Coil F/L	23.1 mm
Coils	9.0
Weight	.95 gm
Maximum Capsule Size	0
Sotax Part Number	6881

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.