


Capsule Sinker Specifications

Sotax (Manual) Style



CAPWST-31

Material	316 Stainless Steel
Finish	Electropolished
Inside Length	31.0 mm
Inside Width	11.0 mm
Wire Diameter	.78 mm
Coil F/L	37.1 mm
Coils	11.9
Weight	3.25 gm
Maximum Capsule Size	000
Sotax Part Number	7349

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.