


Basket Sinker Specifications

Japanese Pharmacopeia Sinker



CUSBSK-JP

Material	316 Stainless Steel
Finish	Electropolish
Inside Length	25 – 26 mm
Inside Diameter	12.0 +/- .2 mm
Wire Diameter	1.0 mm
Pitch	3.0 – 3.5 mm
End Mesh Spacing	3.5 – 4.0 mm
Wire Supports	10
Weight	4.0 gm
Maximum Capsule Size	00

Ordering Information:

QLA Basket Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, or refer to our web site at www.qla-llc.com to see all our product line.