




Low Evaporation Vessel Covers

Economical, Hinged Design



- *Minimizes the evaporation path*
- *Simple low cost configuration*
- *For Distek, Hanson and Vankel baths*
- *Clear Acrylic (consult factory regarding other materials)*
- *Vankel covers utilize the Ease Align system while Hanson and Distek covers register to vessels with locating feet*

Ordering Information:

Catalog Number	Description
COVERD-HG	Cover for Distek Baths
COVERH-HG	Cover for Hanson Baths
COVERV-HG	Cover for VanKel Baths