


Capsule Sinker Specifications

Helix Coil Type



Specifications

	<i>CAPWHT-02</i>	<i>CAPWHT-2S</i>	<i>CAPWHT-04</i>	<i>CAPWHT-4S</i>	<i>CAPWHT-XL</i>	<i>CAPWHT-LG</i>	<i>CAPWHT-XS</i>
Material	Music Wire	316 SS	Music Wire	316 SS	Music Wire	Music Wire	316 SS
Finish	PTFE Coating	Electropolish	PTFE Coating	Electropolish	PTFE Coating	PTFE Coating	Electropolish
Inside Length	.840 in	.840 in	.840 in	.840 in	1.150 in	1.100 in	.885 in
Inside Dia.	.370 in	.370 in	.385 in	.385 in	.465 in	.410 in	.580 in
Wire Dia.	.050 in	.043 in	.043 in	.047 in	.043 in	.043 in	.047 in
Pitch	.156 in	.156 in	.280 in	.280 in	.320 in	.320 in	.156 in
End Coil ID	.175 in	.190 in	.175 in	.190 in	.175 in	.175 in	.190 in
Max Width	.64 in	.64 in	.60 in	.60 in	.69 in	.62 in	.790 in
No. Coils	8.50	8.50	5	5	5.50	6.50	8.5
Weight	2.05 gm	2.85 gm	1.70 gm	1.95 gm	2.15 gm	2.00 gm	3.60 gm
Max Capsule	0	0	0	0	000	000	Special

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.