

IMPURITIES OF DIMETHYLACETAMIDE

Column: **CGB-WAX**, P/N CG-140232

Dimensions: 30m x 0.25mm x 0.25 μ m

Injection: 0.3 μ L, split, 260°C

Carrier gas: H₂, 11 psi (78.8 KPa)

Oven temperature: 75°C(7min) @ 10°C/min to 200°C

Detector: FID, 280°C

Peak Name

- 1- Methanol
- 2- Dimethylformamide
- 3- Dimethylacetamide
- 4- Propylene glycol
- 5- Ethylene glycol
- 6- Monomethyl acetamide

