

columnstore

HPLC column storage

- **Available in 4 sizes**

Columnstore 30, 60, 90 and 120

- **Stackable**

Unique stacking kit allows increase the number of columns you can store at any time.

- **100% Steel Cabinet**

Quality, solid steel construction with chrome-plated D-Ring handles

- **Optional Fittings Tray**

Moulded 6 position tray for easy storage of HPLC tools and fittings

- **Versatile**

Store up to 12 x 12.5cm columns in a single drawer, or a combination of short columns, guard columns and analytical columns. Holds all columns up to 30cm long.



order online @ www.labhut.com

cronus

HPLC columns

● HPLC Columns

Cronus HPLC Columns are manufactured using a variety of different silicas, including Nucleosil, Partisil, LiChrospher, LiChrosorb, Superspher and our own popular Cronusil brands.


All common Phases are available, including C18 (ODS), C8, NH2, Diol and Silica. Special Cronus columns are also available for Peptide and Protein analysis and samples that require an extended pH range.

● Columns & Guard Columns

- Individually Packed
- Individually Tested
- Test Chromatogram Supplied with each Column
- 100% Quality Guaranteed



Authorised Distributor



SMI-Labhut Ltd. Unit D8 Churcham Business Park Churcham Gloucestershire GL2 8AX UK
T : +44 1452 310210 F : +44 1452 300075 sales@labhut.com

All Trademarks Recognised